## English

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**INDICATIONS:** Indicated for Patellofemoral pain syndrome (PFPS) and any condition where additional support to resist lateral patella dislocation is desired. These conditions include Patella Alta, Patella Baja, Patellofemoral Malalignment Syndrome, Patella Subluxation and/or Dislocation, Patellar Tendonitis, Osgood-Schlatter’s Disease, Post-operative Effusion, and Jumper’s Knee.

**WARNING:** For Single Patient Use Only.

**CAUTION:** Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

**CAUTION:** This device will not prevent or reduce all injuries.

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**CAUTION:** Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions. Failure to properly position the brace and fasten each strap will compromise performance and comfort.

**CAUTION:** This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents, or exposure of the hinges to debris or liquids. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

**CAUTION:** Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150°F (65° C.) for any extended period of time.

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

## English

**1. LOOSEN STRAPS**

**(A)** Release the two over-center buckles by pinching the two gray buttons on each of the buckles.

**(B)** Loosen the two remaining straps.

**2. POSITION THE BRACE**

With the patient in a seated position and leg bent at 90°, place the brace on the leg aligning the center of the hinge with the top of the knee cap. Note: Make sure to push the brace posterior so that the hinges run along the mid-line of the leg.

**3. SECURE THE TOP CALF STRAP**

Wrap the top calf strap around the calf as high as possible and loop through the D-ring. Pull the tab until a snug and comfortable fit is achieved and secure.

**4. SECURE THE LOWER THIGH STRAP**

Wrap the lower thigh strap around the thigh and loop through the D-ring. Pull tab until a snug and comfortable fit is achieved and secure.

**5. ADJUST THE CRESCENT FINGER PAD**

Remove the Crescent Finger Pad from underneath the Y-Pad. Relocate the pad to the desired position along the edge of the patient’s knee cap.

**6. SECURE LOWER PATELLO-FORCE STRAP**

**(A)** With the patient’s knee at 20° of flexion, wrap the lower Patello-Force Strap around the leg and secure the over-center buckle.

**(B)** Adjust the strap using the fastener at the strap end until a snug and comfortable fit is achieved.

**7. SECURE UPPER PATELLO-FORCE STRAP**


Repeat Step 6 for the upper Patello-Force Strap.

**8. ADJUST THE PATELLO-FORCE STRAP**

Use the fasteners for additional adjustment.

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